


— REUSABLE PEN · MANUAL PRESS-INJECT

ES03

A reusable, manually driven injection pen — one pen, up to 3 years or 3,000 cycles. Rotate the external dose ring, press to inject, and reuse. Fully customizable, one supply chain: pen, cartridge, needle and carton.

**75 clicks**

Max dose per injection

0.01 mL

Dose step (1 click)

3.0 mL

Cartridge · ISO 11608-2

Manual

Press-inject

3,000 cycles

Reusable · 3 years

1,000 pcs

MOQ

— Specifications

Injection method	Manual — external rotation to set dose, then press the actuator
Dose-setting method	External rotation (turn the external dose ring)
Available doses	60 clicks (0.60 mL) or 75 clicks (0.75 mL)
Dose step	0.01 mL per click
Built-in container	3.0 mL cartridge, ISO 11608-2 compliant
Usage type	Reusable — lifespan of 3 years or 3,000 cycles
Compatible needle	Sterile pen needle, ISO 11608-2 compliant
Body material & colour	Plastic or metal; blue, teal, pink, grey, purple

DIMENSIONS



WHY CUSTOMERS CHOOSE IT

**Reusable by design**

One pen serves for up to 3 years or 3,000 cycles, cutting cost-per-dose and waste.

**Wide cartridge compatibility**

Standard 3.0 mL ISO 11608-2 cartridges across insulin, GLP-1, interleukin and IgE biologics.

— Full customization

**Dose button**

Custom color and molded logo on the button.

**Dose window**

Transparent window printed with your unit of measure (mg, IU or mL).

**Dose scale**

Dial graduations matched to your drug — one click = your set amount.

**Body logo and label**

Screen-printed brand logo and product label on the barrel.

**Cap and body material**

Custom cap and body materials — plastic and metal options.

**Per-part color**

Custom color on clip, body, cap and dial — blue, teal, pink, grey, purple, or your own.

**Hidden needle design**

Optional concealed-needle design for needle-phobic patients.

THE DEVICE

● Five colorways



● Single pen



● Cap · body · cartridge



● Profile detail



Waytop is a wholesale supplier of injection pens and related components. Regulatory market entry, local certifications and compliance are the buyer's responsibility. Waytop does not sell into markets where a given model is subject to an existing authorized distributor or patent. Intended applications, dose ranges and all clinical statements are determined by the final registered labeling for each market.